


More Aesthetics  
More Ethics

# Sterilbeauty

The new autoclave certified for use in beauty centers



Maximum attention  
to the environment



## Small size (7 liters)

The small size allows an easy location within the beauty center rationalizing the spaces in the best way. **The chamber size was designed in order to sterilize a large number of tools and satisfy all the needs.**




## Class S autoclave, for sterilization of wrapped instruments

Instruments must always be packed before sterilization in autoclave, Class S autoclaves are in compliance with directives for packed instruments. **Wrapped instruments are an element of differentiation for the beauty center which wants to communicate hygiene, safety and professionalism.**











## Certified for use in beauty center

This autoclave is certified for use in beauty centers. **Tecno-Gaz is leader in sterilization and pays attention to the world of beauty and aesthetics,** designs and realizes high quality products that can help you in everyday activities.



# Sterilize. Why?

 <h1>1</h1>	 <h1>2</h1>	 <h1>3</h1>	 <h1>4</h1>
<p>► To avoid the risk of infections</p> <p>The parts of the body affected by beauty treatments can host the most harmful of <b>viruses and bacteria that can cause serious disease to humans.</b></p>	<p>► To protect the health of their own collaborators</p> <p>Less disease means more <b>productivity and therefore more earn.</b></p>	<p>► To communicate safety and hygiene to your customers</p> <p>Your customers will ask you <b>wether you sterilize instruments and how.</b></p>	<p>► Only with steam autoclave you can proof a correct sterilization</p> <p>This represents a <b>great legal protection</b> in case of complaints.</p>

 <p>or use in beauty center</p>	 <p>Water saving</p>	 <p>Energy consumption</p>	 <p>Time saving</p>
---	---	---	--

(1 cycle = approx 30 minutes)

# SterilBeauty

Art. AS004ZSM

Class S autoclave



Dimensions (WxHxP)	384x330x530 mm
Dimensions with open door	300 mm
Weight (empty)	27 kg
Weight (full tank and chamber)	31 kg
Acoustic power	< 50 db A
Supply voltage	230 V a.c. +/-10 %
Max power absorbed	1,5 kW
Frequency	50 / 60 Hz
Alimentation cable	2x1,5 mm <sup>2</sup> + hearth

Max working pressure	2.4 bar (relatives)
Max temperature	138 °C
Chamber material	Inox AISI 304
Dimensions (mm)	Ø 164 X 320 mm
Tank volume	2.2 l
Usable tank volume	MAX 1.2÷1.3 l
Cycles	7
Cover material	Polietilene